

Complications cart

CHECKED

Date _____

Initials _____

What to keep on or next to the cart. Check it at every stock check.

PERICARDIUM

- ☐ **Pericardiocentesis tray**
Sheathed needles, J wire, 5 or 6 Fr sheath, multihole pigtail, drainage bottle. Lock pericardiocentesis set and tray, Pericardiocentesis kit
- ☐ **Drainage catheter**
Drainer centesis catheter
- ☐ **Echo machine with a subcostal probe**
In the room for complex PCI and structural cases

FLOW AND THROMBUS

- ☐ **Aspiration thrombectomy catheters**
CAT RX, Export Advance aspiration catheter
- ☐ **Embolic protection**
FilterWire EZ, SpiderFX
- ☐ **IVUS or OCT catheter**
Eagle Eye Platinum, Dragonfly OpStar imaging catheter
- ☐ **Dual lumen microcatheter**
Twin-Pass dual access catheter

SUPPORT, NEARBY

- ☐ **IABP console and catheters**
Cardiosave IABP (Hybrid and Rescue), Sensation Plus IAB catheter
- ☐ **Impella CP and its sheath**
Impella CP with SmartAssist
- ☐ **Venoarterial ECMO circuit and cannulas**
Cardiohelp system
- ☐ **Distal perfusion sheath**
Super Arrow-Flex sheath

PERFORATION

- ☐ **Covered coronary stents in the sizes**
PK Papyrus 2.5 to 5.0 mm; Graftmaster 2.8 to 4.8 mm. PK Papyrus, Graftmaster RX
- ☐ **Coils, 0.014 and 0.018 inch compatible**
Not in the directory yet
- ☐ **Coil compatible microcatheter**
Finecross M3
- ☐ **Semicompliant balloons to seal, vessel sized**
Trek and Mini Trek
- ☐ **Perfusion balloon**
Ringer perfusion balloon catheter
- ☐ **8 Fr guide catheters**
Launcher guide catheter

VASCULAR

- ☐ **Large bore sheaths**
Sentrant introducer sheath, DrySeal Flex introducer sheath
- ☐ **Peripheral balloons to tamponade a bleed**
Armada 35, Mustang
- ☐ **Aortic occlusion balloon**
Coda balloon catheter
- ☐ **Peripheral covered stents**
Viabahn VBX, Viabahn
- ☐ **Crossover sheath**
Pinnacle Destination guiding sheath
- ☐ **Vascular plug for a bleeding branch**
MVP micro vascular plug
- ☐ **Stiff and hydrophilic 0.035 inch wires**
Amplatz Super Stiff, Glidewire
- ☐ **Closure devices and compression**
Perclose ProStyle, MANTA vascular closure device, FemoStop Gold, TR Band

DRUGS

- ☐ **Atropine**
- ☐ **Epinephrine, 1 mg/mL and 0.1 mg/mL**
- ☐ **Vasopressin, norepinephrine, phenylephrine**
- ☐ **Protamine**
Dose: Cath Med App
- ☐ **Amiodarone and lidocaine**
Lidocaine hydrochloride injection
- ☐ **Adenosine**
Dose: Cath Med App
- ☐ **Nicardipine, nitroprusside, verapamil**
Verapamil hydrochloride injection
- ☐ **Nitroglycerin**
Nitroglycerin in dextrose
- ☐ **Cangrelor, eptifibatide or tirofiban**
Kengreal (cangrelor), Aggrastat (tirofiban). Dose: Cath Med App
- ☐ **Diphenhydramine, H2 blocker, hydrocortisone, albuterol, glucagon**
- ☐ **Naloxone and flumazenil**
- ☐ **Dobutamine**
Dose: Cath Med App
- ☐ **Thrombin for pseudoaneurysm**

RETRIEVAL

- ☐ **Coronary snares, 2 to 7 mm loops**
Not in the directory yet
- ☐ **Large multi loop snares, 18 to 30 mm**
Not in the directory yet
- ☐ **Retrieval forceps or bioptome**
Not in the directory yet
- ☐ **Small balloons, 1.0 to 1.5 mm**
Sapphire II Pro 1.0
- ☐ **Guide extension and trapping device**
GuideLiner V3, TrapLiner
- ☐ **Long sheaths for device retrieval, 12 to 16 Fr**
DrySeal Flex introducer sheath

RHYTHM

- ☐ **Defibrillator with pads**
- ☐ **Temporary pacemaker generator**
Medtronic single chamber temporary pacemaker
- ☐ **Pacing catheter kit and a venous sheath**
Arrow pacing catheters and kits
- ☐ **Pacing guidewire**
Wattson temporary pacing guidewire
- ☐ **Mechanical CPR device**